

## Sample Diluent Buffer (1X – Ready to use)

Sample diluent that minimizes non-specific binding, cross-reactivity and matrix interference.

Reference	IMS0514
Volume	10 ml
Storage:	2-8°C
pH-value:	7.2 ± 0.2
Preservative:	contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)
Expiry date when stored unopened:	please refer to the label on the bottle

**\*For research use only, not for diagnostic use**

### Instructions for use

- Before use, temper the buffer between +18°C y +24°C (room temperature) for 30 minutes.
- Samples are diluted 1:20 sample diluent buffer (1X) and mixed with a vortex agitator according the below figure. Diluted samples should be assayed within 8 hours.

