

## Wash Buffer (10X)

PBS 10% BSA.

Reference:	IMS0515
Volume:	20 ml
Storage:	2-8°C
pH-value:	7.4 ± 0.2
Preservative:	contains 10 % albumin in 10mM sodium phosphate, 150 mM NaCl, pH 7.4, contains KATHON anti-microbial agent.
Expiry date when stored unopened:	please refer to the label on the bottle

**\*For research use only, not for diagnostic use**

## Instructions for use

- Before use, temper the buffer between +18°C y +24°C (room temperature) for 30 minutes.
- If crystallization is observed in the concentrated buffer during storage, warm to 37°C and properly shake before diluting.
- Dilute the contents of the 10X assay buffer to 1X (PBS) in PBS 1X, for use flow cytometry assays (Exosomes and SARS-CoV-2 Immunostep Assays). To carry out the dilution, the amount required for the assays is removed from the concentrate bottle and diluted 1:10 in PBS 1X.

